

INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No.	Report No.:	Page No.
		1 of 21

SERIAL NO.	ITEM DESCRIPTION	PAGE NO.
1.0	PROTOCOL APPROVAL	2
2.0	OVERVIEW:	3
2.1	Objective	3
2.2	Purpose	3
2.3	Scope	3
2.4	Responsibility	4
2.5	Execution Team	4
3.0	ACCEPTANCE CRITERIA	5
4.0	REVALIDATION	5
5.0	INSTALLATION QUALIFICATION PROCEDURE	6
5.1	System Description	6-8
5.2	Instruction for Filling the Checklist	9
5.3	Installation Check-List	10
5.4	Identification of Major Components	11-12
5.5	Verification of material of construction	13
5.6	Identification of Supporting Utilities.	14
5.7	Identification of Safety Feature(s)	14
5.8	Identification of Standard Operating Procedure	15
5.9	Test instrument details	15
5.10	Identification of Component to be Calibrated	15
5.11	Verification of Drawing and Documents	16
5.12	Abbreviations	17
5.13	Deficiency And Corrective Action(s) Report(s)	18
5.14	Annexure(s)	19
6.0	INSTALLATION QUALIFICATION FINAL REPORT	20
6.1	Summary	20
6.2	Conclusion	20
6.3	Final report approval	21



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		2 of 21

1.0 PROTOCOL APPROVAL:

Signing of this approval page of Protocol indicates agreement with the qualification approach described in this document. If modification to the qualification approach becomes necessary, an addendum shall be prepared and approved. The protocol cannot be used for execution unless approved by the following authorities.

This Installation Qualification protocol of Auto Coater has been reviewed and approved by the following persons:

FUNCTION	NAME	DEPARTMENT	SIGNATURE	DAT E
PREPARED BY		QUALITY ASSURANCE		
REVIEWED BY		PROJECTS / ENGINEERING		
REVIEWED BY		PRODUCTION		
APPROVED BY		QUALITY ASSURANCE		



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		3 of 21

2.0	OVERVIEW:
0.4	OD IEOTIVE
2.1	OBJECTIVE:
	The objective of developing and executing this protocol is to collect sufficient data
	pertaining to the Auto Coater and define the installation qualification requirements and
	acceptance criteria for the Auto Coater. Successful completion of these installation
	qualification requirements will provide assurance that the Auto Coater was installed as
	required in the manufacturing area.
2.2	PURPOSE:
	The purpose of this protocol is to establish documentary evidence to ensure that the Auto
	Coater system received matches the Design specification and also to ensure that it is
	properly and safely installed.
2.3	SCOPE:
	This Protocol is applicable to installation of Auto Coater at the tablet manufacturing facility
	at (Company Name) & the subsequent documentation.
2.4	RESPONSIBILITY:
	The following shall be responsible:
	Quality Assurance officer/ Executive - Preparation of protocol
	Projects / Engineering Head – For execution
	Production Head – For execution support

Quality Assurance Head – For adequacy and final approval

2.5 **EXECUTION TEAM:**



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:		
Protocol Reference No.	Report No.:	Page No.
		4 of 21

The satisfactory installation of the Sejong Auto Coater shall be verified by executing the qualification studies described in this protocol. The successfully executed protocol documents that the Auto Coater is installed satisfactorily.

Execution team is responsible for the execution of installation of Auto Coater Execution team comprises of:

DEPARTMENT	DESIGNATION	NAME	SIGNATURE	DATE
PROJECTS/ ENGINEERING				
PRODUCTION				
QUALITY ASSURANCE				



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	ter Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		5 of 21

3.0	ACCEPTANCE CRITERIA:
3.1	The Auto Coater shall meet the system description given in design qualification.
3.2	The Auto Coater shall meet with the acceptance criteria mentioned under the topic
	"Identification of major components"
3.3	The Auto Coater system shall be operated by PLC.

4.0	REVALIDATION CRITERIA:	
	The Auto Coater has to be revalidated if	
	There are any major changes in system components which affect the	
	performance of the system	
	After major breakdown maintenance is carried out.	
	As per revalidation date and schedule	



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.

-	_
 	6 of 21

5.0	INSTALLATION QUALIFICATION PROCEDURE:

5.1 AUTO COATER (Modal)SYSTEM DESCRIPTION:

Process Equipment Description

The purpose of Sejong Auto coating Machine SFC 100FN is to coat tablets of all common sizes and shapes whether it is for film or aqueous coating. The hot air in this equipment will enable a very efficient heat and mass transfer. The hot air flow is controlled by PLC and is programmable to accommodate for various types of tablets and their coating materials.

Complete machine can be divided in following sub sections:

- Inlet System
- Pan
- Spray System
- PLC System
- Cleaning system
- Exhaust System
- Control System

Inlet System: It consists of pre & medium air filters, HEPA filters (Optional), air heating unit and airflow & temperature measuring devices. Inlet air handling unit (Inverter control) receives air from the environment and processes the inlet air and controls the inlet airflow. **Pan:** It consists of the semi-perforated fully enclosed pan which contains 8 ducts that engages with the incoming air and the outgoing air through arc shaped inlet duct from the heat exchanger and outlet duct that goes to the exhaust. At any point of time two circular cross sectional area of the pan engages with the hot air inlet and outlet.



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		7 of 21

Spray System. It consists of peristaltic pump in case of film coating. This pump is selected based on the viscosity of the solution that has to be coated. This pump is eventually connected to the tubing that leads to high-pressure low volume guns. The angle of spray and the flow rate are externally controllable.

A programmable logic controller (PLC) will execute the direction of the sequence of operation. The PLC of OMRON (or equivalent) is used for this purpose.

5.2 PLC Description

The main function of a PLC is to translate the instructions into the digital or analog codes needed to operate the device or machine.

PLC collects data from field instrumentation & display the information on the operator station. The instruments are connected to the system equipment and piping. The collected data will be utilised by the PLC for process control.

The user interface, based in an industrial type HMI, will assist the operator to supervise and control the process. Based on the displayed information the operator, by means of the user interface, can provide commands to the PLC.

The PLC then executes the operator instructions. An OMRON PLC CQM1H-CPU44 has been chosen as the Central Processing Unit (CPU). There are six inner boards available, which are mounted in Slot 1 or Slot 2 of a COMQ1H-CPU44. The protocol macros are made on the CX – Protocol Support Software and then recorded in the Serial Communication Board, where they can be executed any time using the PMCR instruction in the CPU unit's ladder diagram.


The PLC system layout for the coating machine automation is as shown below:



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No.	Report No.:	Page No.	
		8 of 21	



Operating Terminal

HP Laserjet/Inkjet Printer



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		9 of 21

5.2	INSTRUCTION FOR FILLING THE CHECKLIST
5.2.1	In case of the compliance of the test use the word 'Complies' otherwise use 'Does not
	comply' to indicate non-compliance.
5.2.2	For identification of the components of the equipment and utilities use the word "yes" to
	show its presence and use 'No' to indicate the absence of the identity
5.2.3	Give the detailed information in the summary and conclusion part of the Installation
	Qualification report.
5.2.4	Whichever column is blank or not used 'NA' shall be used.



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		10 of 21

5.3	INSTALLATION CHECKLIST:
	Installation checklist is as follows:

Sr.NO.	STATEMENT	YES / NO	CHECKED BY (SIGN)	DATE
1.	Verify that the "As Built" drawing is complete and represents the design concept.			
2.	Verify that major components are securely anchored and shock proof.			
3.	Verify that there is no observable physical damage.			
4.	Verify that there is sufficient room provided for servicing.			
5.	Verify that all piping and electrical connections are done according to the drawings.			
6.	All access ports are examined and cleared of any debris.			
7.	Safe electrical connections.			
8.	Sufficient room provided for maintenance.			
9.	Equipment identification nameplate visible.			
10.	Units installed on foundation are secure in place as per manufacturer's recommendations.			



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:		
Protocol Reference No.	Report No.:	Page No.	
		11 of 21	

5.4 IDENTIFICATION OF MAJOR COMPONENTS:

Describe each critical component and check them and fill the inspection checklist.

System Components	Design Sp	ecification	Complies / Does Not Comply	Checked By (Sign)	Date
	Model	Modal Name			
1. Equipment	Output Capacity	120 - 150 Kg			
Description	Power of Main motor	3.7KW			
	Hot water Pressure	3 Kg/cm sq.			
	Volume of Heating	42000 K Cal			
	Pan diameter	1300mm			
	Material	SS 316L			
2. Pan	Spray Sliding Cylinder	SS 304			
	Rabbit Ear Baffles	7 Nos.			
	Pan Motor	3 HP			
	Pan speed	1 - 12 RPM			
3. Solution	Capacity	100 L			
Tank	Material of construction	SS316L			
	Manufacturing Design	Casters mounted.			
4. Air Device	Main Air Regulator	0~10 Kg/cm sq.			
Rating	Atom Air Regulator	0~10 Kg/cm sq.			
	Spray Gun Air Regulator	0~10 Kg/cm sq.			
	Air Pump Regulator	0~10 Kg/cm sq.			
	Main Pressure Gauge	0~10 Kg/cm sq.			
	Air Flow Meter	0~500 Lit/min.			
	Manometer	0~500 mm Aq			
5. Pump	Manufacturer	XYZ			



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No. Report No.: Page No.

--- 12 of 21

System Components	Design Specification		Complies / Does Not Comply	Checked By (Sign)	Date
	Model	504S/RL			
	RPM	Max220			
	Shaft Torque	2.2 Nm			
6. Other	Touch Screen	XYZ			
Components	Proximity switches	Autonics Corporation			
	Flowmeter	Dwyer			
	Magnehelic Differential pressure gauge	Dwyer			
	Solenoid Valve	TPC Pneumatics			
	VFD	Mitsubishi Electric			
7. Spray Gun	Manufacturer	XYZ			
	Model	1AV-HVLP			
8. Exhaust	Manufacturer	XYZ			
Blower Air Fan Motor	Speed Controller	XYZ VFC			
INIOIOI	Filter	Polyester filter bag with auto shaking function			
9. Operating Control Panel	PLC Make	XYZ			
	PLC Model	Sysmac CJ1G			
	PLC CPU	CPU 44			



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No. Report No.: Page No.

--- 13 of 21

5.5 Verification of Material of construction:

Sr. No.	Name of components	Material of construction	Method of verification	Verified by (Sign & Date)
1	Rotary Inlet Outlet Valve	SS 41		
2	Outside Door	SS304		
3	Outside Cover	SS304		
4	Pan Door	SS316 & Strengthened Glass		
5	Auto damper	SS304		
6	Coating Pan	SS316		
7	Film Tank	SS 316		
8	Coating Door	SS 304		
9	Machine Body	Frame: SS41, Cover: SS 304		
10	Liquid Tank Flange	SS304		
11	Tank Top Cover	SS304		
12	Liquid Tank	SS316		
13	Liquid Tank Oil Seal Cover	SS304		
14	Liquid Tank Impeller shaft coupling	SS304		

5.6 IDENTIFICATION OF SUPPORTING UTILITIES:



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		14 of 21

UTILITY	PROPERLY IDENTIFIED & CONNECTED (YES/NO)	CHECKED BY (SIGN)	DATE
1) Electricity:			
3 Phase, 415Volts, 50Hz			
2) Compressed air			
6 kg/cm ²			
12 mm, 500 L/min			
3) Hot water			
65000 Kcal/Hr, 2Kg/cm2			

5.7	IDENTIFICATION OF SAFETY FEATURES:			
	Identify and record the safety features (if any) and their function in following tables:			wing tables:
Safety Features Description		Function	Identified By (Sign)	Date
Earthing	1	To avoid electrical shocks due to		
		leakage current.		
Alarm Message				
1. Purging Air Pressure		If the air pressure drops below the		
Low		minimum set level, then this alarm		
		will trip the machine.		
2. Machine Stoppage		If emergency stop is operated the		
		alarm will generate		
3. Main Motor Overload		If the Sharft motor is overload the		
		alarm will generate and trip the		
		process		

5.8 IDENTIFICATION OF STANDARD OPERATING PROCEDURE (SOP)



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No.	Report No.:	Page No.
		15 of 21

5.9	TEST INSTRUMENT DETAILS
	Instrument/Equipment required:
	1.Tacho Meter
	2. Anemometer
	3.Temparature sensor
	4.Pressure Gauge
	5. Photometer

5.10 IDENTIFICATION OF COMPONENT TO BE CALIBRATED

In the Auto Coater, following are the components, which needs calibration. Following are the components:

- 1.Pressure Measurement (Purging Air)
- 2. Temperature Measurement
- 3. R.P.M.

5.11	VERIFICATION OF DRAWING AND DOCUMENTS:
	Following documents are reviewed and attached as listed below:



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No.	Report No.:	Page No.
		16 of 21

Sr. No.	DRAWING AND DOCUMENT DETAIL	CHECKED BY (SIGN)	DATE



INSTALLATION QUALIFICATION

Equipment Name: Auto coater Equipment No.:

Protocol Reference No. Report No.: Page No.

--- 17 of 21

5.12	ABBREVIATIONS
	Following Abbreviations are used in the installation qualification protocol of Air handling
	unit
	MOC: Material of construction
	RPM: Rotation per minute
	Nm: Neuton meter
	PLC: Programming Logic Controller
	ACT : Auto Coater

DEFICIENCY AND CORRECTIVE ACTION (S) REPORT (S)



INSTALLATION QUALIFICATION

INSTALLATION QUALIFICATION		
Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		18 of 21

Following deficiency was identified and corrective actions taken in consultation with the validation
team.
team.
Description of deficiency:
Corrective action(s) taken:

Reviewed By:

Date

5.13

5.14 Annexure(s):



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.
		19 of 21

Sr.No.	Annexure No.	Title of Annexure

PHARMAGUDDU.COM INSTALLATION QUALIFICATION Equipment Name: Auto coater Protocol Reference No. Report No.: Page No. --- 20 of 21

All the IQ data sheets and discrepancy report shall be reviewed by validation team to prepare summary report. The summary of IQ shall be used to draw conclusion for approval of installation qualification report.

6.1	SUMMARY
6.2	CONCLUSION

6.3 FINAL REPORT APPROVAL



INSTALLATION QUALIFICATION

Equipment Name: Auto coater	Equipment No.:	Equipment No.:	
Protocol Reference No.	Report No.:	Page No.	
		21 of 21	

It has been verified that all tests required by this protocol are completed, reconciled and attached to this protocol or included in the qualification summary report. Verified that all amendments and discrepancies are documented, approved and attached to this protocol.

Signature in the block below indicate that all items in this qualification report of Auto Coater have been reviewed and found to be acceptable and that all variations or discrepancies have been satisfactorily resolved.

NAME	DESIGNATION	DEPARTMENT	SIGNATURE	DATE
		PROJECTS / ENGINEERING		
		PRODUCTION		
		QUALITY ASSURANCE		